

A novel and rapid process for ambient temperature preservation of microorganisms

Sebastien Keller - PhD – Technological Development Manager – <u>keller@imagene.eu</u> - www.imagene.eu

Imagene is a french biotech company

- Imagene has developped a new technology for preservation of biological samples at **ambient temperature.**
 - Dedicated to long term storage (years), and in particular to Biobanking
 - Offers an alternative to the constraints of cold storage

The technology is already validated for DNA and RNA



What is behind Imagene technology?

Imagene has demonstrated with DNA that storage at room temperature is feasible over long period of time if the samples are protected from degradation factors, mainly Oxygen and Moisture:



Redrawn from Colotte et al. Biopreservation and Biobanking 2011

The imagene process

imagene

Imagene solution : Storage of biospecimens in a patented container: paperclip-size, airtight **minicapsules**, away from all degradation factors.



Main Advantages of Imagene Technology

Long term stability demonstrated for DNA and RNA 1 cut / 10 000 nucleotides / century for DNA and 1 cut / 700 nucleotides / century for RNA

(Estimated from accelerated aging studies)

Compared to cold storage:

- Substantial reduction of required storage space (up to 200 000 minicapsules in 6 m3)
- Substantial reduction of storage cost (no electric consumption, no maintenance...)
- Compared to freeze-drying
 - Faster process and easier to automate
 - Metallic container : no risk of breakage during transport



The Anvbis³ project

Project started in 2013.

Objectives:

- Evaluate the potential of the Imagene's process to be used for room temperature long term preservation of bacterial strains.
- Aim is to preserve enough viable bacteria to allow culture growth after several years of room temperature storage.

Methodology

- Use of 2 model bacteria: *Pantoea dispersa* (freeze drying resistant)
 & *Aeromonas salmonicida* (freeze drying sensitive)
 - Performance criteria : Viability Rate VR expressed as:

 $VR(\%) = \frac{Log(viable concentration after processing)}{Log(viable concentration before processing)} x \ 100$



The imagene process

imagene

Imagene solution : Storage of biospecimens in a patented container: paperclip-size, airtight **minicapsules**, away from all degradation factors.



The Anvbis³ project – drying step

- Different Vacuum drying conditions were tested (product formulation, drying parameters...)
 - Best conditions for the two strains were compared with a standard freeze-drying protocol, performed both in glass vials and in minicapsules :



The Anvbis³ project – encapsulation step

Encapsulation does not significantly alter the viability rate.





The Anvbis³ project – storage step

Study of storage at room temperature and 37°C in minicapsules



imagene

The study is on-going (RT).

Accelerated aging modelling is an issue

The Anvbis³ project – preliminary conclusion

- Imagene's technology is compatible with the storage of bacterial strains.
- More optimizations and characterizations are needed
 - Long term measurement
 - Test of other strains
 - Control of other strain functions
 - Aging modelling

....

In parallel, development of a prototype platform for drying and encapsulation of class II & III bacterial and viral strains.



Anvbis³ encapsulation platform

Platform for safe encapsulation of class II and III microorganisms.

image

- Include all the steps of the Imagene Process
 - Semi automated
 - Up to 100 samples per day





Acknowledgements:

- DGA (French Defense procurement agency – funder of the project)

- Pasteur Institute (partner)
- Genopole Evry





THANK YOU



keller@imagene.eu

13

